

SVCell technical overview

Automated Image Recognition with Recipes

- Single click application execution provides automated image import, enhancement, segmentation, measurement and classification
- Individual recipes can be ordered in a recipe list for each application and saved
- Easy recipe sharing and porting
- Application images can be loaded into a Field of View (FOV) list for batch processing with recipes
- The entire experimental outcome can be quickly evaluated inside the unique RFrame integrated image and data analysis user interface

Image Viewing

- 4 dimensional (x, y, channel, time) FOV viewing and processing
- N-channel color viewing and processing
- Time based FOV navigation
- FOV thumbnail view
- FOV merging by channel or time
- Zooming and magnification tools
- Confidence map viewing panel
- Unlimited number of masks for each channel
- Intensity mapping for visualization
- Time-lapse image alignment

Data Analysis

- Large number of object and field measurements
- Spatial and intensity calibrations
- Object statistics displayed in spreadsheet
- Region of interest (ROI) graphics overlain on the objects in images to facilitate selection linking
- Selection linking between object spreadsheet rows, image ROI overlays and chart data points to facilitate integrated image and data review
- Dot plot and histogram charts
- Trace-plots
- Measurements, statistics and field measurements displayed in the RFrame spreadsheet and charts
- Gate from charts to create sub-populations
- Display sub-population data points overlain on the parent population charts
- Color ROIs and data by FOV and sub-population
- Spreadsheet sorting
- Export data to comma separated values format (MS Excel readable)

Import Recipe Teaching

- Automate the import of individual frames on the filesystem into a multi-channel, time-lapse FOV list
- Relative import option to support application deployment

Segmentation Recipe Teaching

- Unique "teach by example" drawing interfaces
- Two modes of confidence mapping
- Confidence map viewing
- Threshold tools with bimodal offset option, as well as histogram percentile units option
- Homotopic threshold provides automatic partition of adjacent objects
- Adaptive threshold provides automatic image normalization prior to thresholding
- Fill holes & remove objects operation
- Overlap separation operation
- Guided partition operation
- User teaching, threshold and operation options are automatically saved to a segmentation recipe
- Unlimited number channels and masks can be processed with a single segmentation recipe
- Recipe can be easily updated through additional teaching; update teaching is fast
- Recipe performance is stable and predictable

Enhancement Recipe Teaching

- Automate common tasks such as renaming masks and channels, deleting masks and channels, image bit depth conversion, image copying, etc.
- A large set of image processing functions covering arithmetic, bitwise, comparison, distance transforms, common filters, geometric, logical, and zone of influence operations
- An advanced morphological processing library containing twenty-six morphological filters
- Grayscale and binary labeling functions for removing objects based on metrics such as size, mean intensity, max intensity, standard deviation and compactness
- Recipe can be easily created by selecting and saving operations from the processing history

Measurement Recipe Teaching

- Allows creation of large number of field and object measurements
- Inter channel measurements
- Support for inter-channel, multi-mask objects

SVCell technical overview

- Ratio measurements
- Cell templates for easy measurement

Decision recipe teaching

- Easy to use, teach by example labeling interface
- Sophisticated classification of objects using SVision regulated decision tree technology
- Automatic creation of the decision recipe containing the decision tree and user teaching
- Transparent tree shows all the decision rules
- Decision tree contingency table
- Decision tree parallel coordinates plot
- Object linking supported inside the recipe between the tree, images and charts
- Recipe classifies objects which are then displayed as subsets in the RFrame spreadsheet
- Classified objects can be saved as a labeled mask
- Recipe can be easily updated through additional teaching; update teaching is fast
- Recipe performance is stable and predictable

Memory Management

- User specifiable memory management, large image mode, and image database performance configurations

Operating System and File Format Support

- Windows® XP, XP Pro and Windows® Vista
- Read and write BMP, DIB, GIF, JPG, JPEG, TIF, TIFF, PNG, RAW, RLE
- Read and write multi-dimensional formats such as OME (Open Microscopy Environment), ICS / IDS, STK and multi-page TIF
- FOVs can be saved to OME compliant SVCell TIF, which saves the FOV dimensionality along with masks, object ROIs, data and inherent relations
- 8, 16, and 24-bit grayscale and color images are supported

Recommended Platform

- Processor: Intel Pentium 4 or later
- Processor speed: 3.0+ GHz
- RAM: 2+ GB high speed memory
- Operating System: Windows XP, Service Pack 2

	<i>Fixed Point configuration</i>	<i>Live Cell configuration</i>
All other SVCell features	✓	✓
Load multi-channel FOV	✓	✓
Load time-lapse FOV		✓
Create, import, enhancement, measurement and segmentation recipe	✓	✓
Create decision recipe	Option	Option
Create tracking recipe	Option	Option
Load recipe	Option	Option
Save recipe	Option	Option



3633 136th PI SE Suite 300
Bellevue, WA 98006, USA
www.svisionllc.com
www.svcell.com

SVision LLC is an emerging leader in the application of learning technologies to practical imaging applications. SVision is an innovator with 30 issued and 18 pending U.S. patents. SVision software is deployed in advanced semiconductor wafer inspection systems that can be taught for new applications without technical intervention. SVision staff have over 140 years of practical imaging solution experience, 152 technical publications, 82 issued U.S. patents, and a history of technical breakthroughs.

SVCell is protected, without limitation, by U.S. Patents 6400849, 6404934, 6456741, 6463175, 6504959, 6507675, 6640008, 6859550, 6941288, 7031529, 7031948, 7096207, 7110603, 7133560, 7139764, 7142718, 7149357, 7203360, 7233931, 7263509, 7293000